

Device Schedule:

Intended Purpose as per the Instructions for Use:

VICRYL RAPIDE™ Suture is intended for use in soft tissue approximation where only short term wound support is required and where the rapid absorption of the suture would be beneficial. Due to its absorption profile VICRYL RAPIDE™ Suture is useful for skin closure, particularly in paediatric surgery, episiotomies, circumcision and closure of oral mucosa. VICRYL RAPIDE™ Suture is also successfully used in ophthalmic surgery for conjunctival sutures.

Device Name	Model	Type (Codes as per (EU) 2017/2185)	Risk Classification	Basic UDI-DI
VICRYL RAPIDE (Polyglactin 910) Synthetic Absorbable Suture	Vicryl Rapide	MDN 1104	Class III, Implantable	0705031a012094L

Additional Information:

Surgical needle and suture combinations from within the following limits define the device range for VICRYL RAPIDE™ (Polyglactin 910) Synthetic Absorbable Suture

Suture Characteristics	Range
Suture Material (Absorbable/Non-Absorbable)	Absorbable
Suture Gauge Size	0.4 – 4.0 (Metric)
Suture Length	45 cm – 150 cm
Suture Dyed/Undyed	Dyed / Undyed
Suture Color (If dyed)	Violet
Coated/Uncoated	Coated
Multifilament/Monofilament	Multifilament
Accessories to suture type	N/A
Needled/Non-Needled	Needled / Non-Needled
Number of Needles per Suture	Single Armed / Double Armed
Needle Material	420 SS, 4310 SS, 4310 FB SS and ETHALLOY
Needle Coating	Silicone, MULTIPASS, CERBERUS
Needle Shape	Straight / Curve
Needle Color	Silver / Black
Needle Length	5.5 mm – 60 mm
Needle Wire Diameter	0.15 mm – 1.27 mm