

Device Schedule:

Intended Purpose as per the Instructions for Use:

LAPRA-TY™ II Clips were developed in order to secure the suture ends of a coated VICRYL™ / VICRYL™ Plus Suture of size USP 2/0, 3/0, or 4/0 (metric sizes 3, 2, or 1.5) after application of a tension-free suture. This ensures a reliable approximation for about 14 days.

Product Name: LAPRA-TY™ II Clips

Basic UDI-DI Device Name: LAPRA-TY II Suture Clips

Basic UDI-DI Device Model: LAPRA-TY II Suture Clips

Basic UDI-DI: 0705031a000083Z

Risk Classification: Class III Implantable

Type (Codes as per (EU) 2017/2185): MDN 1104

Product Code	Product Description	Number of Clips Per Cartridge	Number of Cartridges Per Box
MIC4032	3mm wide LAPRA-TY II Clips in Cartridge	2	6
MIC4034	3mm wide LAPRA-TY II Clips in Cartridge	4	6
MIC4036	3mm wide LAPRA-TY II Clips in Cartridge	6	6